



Estd. 1964

YSI-106

Automatic Slide Staining Machine



Suitable to work on 220V, single phase, 50 Hz, AC supply.

12 Station

Manual Slide Staining procedure occupies considerable amount of technician's time. To accelerate the process "Unimeditrek Slide Staining Machine" was designed on modern principles for carrying out practically all slide Staining techniques of Cytology & Histology.

Unimeditrek Slide Staining Machine

"12 Stations", is compact and rugged equipment designed on latest technology for automatic staining of slides in different reagents according to preprogrammed chronological user selectable cycles. Throughout the process the slides are shaken (Switch selectable) for effective staining in a slide (Capacity 16 Slide). The last trough i.e. 12th Station is provided with an inlet and outlet to enable washing of the residual chemicals.

The unit is supplied complete with 11 Glass Staining troughs, 1 trough with inlet and outlet, 1 Slide carrier, 5 Timing Discs for chronological programming calibrated for 60 minute and 1 notch cutting plier.



Manufactured By:

Manual Stainer Slide Staining Set

Simple and cheap system to use for cytological and histological staining. The set includes a twelve places stainless steel holder for dishes, twelve solvent resistant dishes with lids and a slide a three place stainless steel holder, three dishes and a slide staining holder.

	Description	Dishes
YSI-109(HS)	Hoematological set	3
YSI-109	Slide Staining set	12



Accessories and Refill

Description	Dishes
Stainless Steel holder - 12	1
PEPT dish with lid	12
TPX transparent dish with lid	12
Slide staining holder for 25 slides	6
Slide staining holder with plastic handle	6
Double slide staining holder (for 50 slides)	1

